The background of the slide is a soft-focus photograph of medical supplies. On the left, a white plastic pill bottle is tipped over, with several light blue, oval-shaped tablets scattered on the surface. To the right, a silver stethoscope with a black tube is resting on a white computer keyboard. The overall lighting is bright and clinical, with a shallow depth of field that keeps the text area sharp while blurring the background elements.

Brain Tumor Classification in MRI Images

Derya Kurt
June 2021

01

INTRODUCTION

05

MODEL EVALUATION

02

PROBLEM STATEMENT

06

CONCLUSIONS

03

DATA AND PREPROCESSING

07

RECOMMENDATION

04

CLASSIFICATION METHODS AND MODELS

08


FUTURE WORK



01

INTRODUCTION

BRAIN TUMOR





02

PROBLEM STATEMENT

Magnetic Resonance Imaging (MRI)

Most common technique for tumor
diagnosis

Limitations of MRI

- Time Consuming during MRI Analysis
- Accuracy depends on radiologists experience

Objective

Develop a Deep Learning approach to detect and
classify Brain Tumor in a short time





03

DATA AND PREPROCESSING

DATA OVERVIEW

2,870

MENINGIOMA

822

GLIOMA

826

PITUITARY

827

NO TUMOR

395

DATA PREPROCESSING

IMAGE PROCESSING

SIZE Different Image Shapes (250 X 250) (130 X 130) (512 X 512) (300 X 236)
Images Resized to 150 X 150

COLOR RGB ---> GRAYSCALE

PIXEL VALUE Normalized to [0,1]

CLASS LABELS ['glioma', meningioma, no tumor, pituitary'] $\xrightarrow{\text{LabelEncoder()}}$ [0, 1, 2, 3]


OneHotEncoding

Glioma	[1 0 0 0]
Meningioma	[0 1 0 0]
No Tumor	[0 0 1 0]
Pituitary	[0 0 0 1]


DATA PREPROCESSING

DATA AUGMENTATION

- Avoids Overfitting
- Shear - Zoom - Width & Height Shift - Horizontal - Flip Rotation



Original Image



Augmented Image



04

CLASSIFICATION METHODS AND MODELS

Multi-Class Classifiers Method

CNN Model Architecture

Layer (type)	Output Shape	Param #
conv2d_10 (Conv2D)	(None, 150, 150, 64)	640
max_pooling2d_10 (MaxPooling)	(None, 75, 75, 64)	0
conv2d_11 (Conv2D)	(None, 75, 75, 128)	73856
max_pooling2d_11 (MaxPooling)	(None, 38, 38, 128)	0
conv2d_12 (Conv2D)	(None, 38, 38, 128)	147584
max_pooling2d_12 (MaxPooling)	(None, 19, 19, 128)	0
conv2d_13 (Conv2D)	(None, 19, 19, 128)	147584
max_pooling2d_13 (MaxPooling)	(None, 10, 10, 128)	0
conv2d_14 (Conv2D)	(None, 10, 10, 128)	147584
max_pooling2d_14 (MaxPooling)	(None, 5, 5, 128)	0
flatten_2 (Flatten)	(None, 3200)	0
dense_4 (Dense)	(None, 512)	1638912
dense_5 (Dense)	(None, 4)	2052
=====		
Total params: 2,158,212		
Trainable params: 2,158,212		
Non-trainable params: 0		

5 – Convolutional Layer

5 – MaxPooling Layer

1 - Fully Connected Layer

Output Layer with 4 classes

.....



Multi-Class Classifiers Method

Model Compiling		Model Training CALLBACKS
OPTIMIZER	ADAM	EARLYSTOPPING
LEARNING RATE	0.01	Stop training when the monitored metric validation set accuracy stops improving
LOSS FUNCTION	CATEGORICAL CROSSENTROPY	MODELCHECKPOINT
METRICS	ACCURACY	Save the model weights at the maximum of the monitored quantity

Multi-Class Classifiers Method


CNN Multiclass Model Evaluation

Validation Accuracy - **0.8803** Test Accuracy - **0.8815**



Binary-Classifiers Method

Meningioma Binary Classifier



Model Architecture

Same architecture with Multiclass CNN Model:

- 5 – Convolutional Layer
- 5 – MaxPooling Layer
- 1 - Fully Connected Layer
- Output Layer with 2 classes

Model Compiling



Same paramaters with Multiclass CNN Model

Model Training

Binary-Classifiers Method


Meningioma Binary Model Evaluation

Validation Accuracy - **0.8243** Test Accuracy - **0.8061**



Binary-Classifiers Method

Glioma Binary Classifier



Model Architecture

3 – Convolutional Layer
3 – MaxPooling Layer
1 - Fully Connected Layer
Output Layer with 2 classes

Model Compiling

Model Training



Same paramaters with Multiclass CNN Model

.....

Binary-Classifiers Method


Glioma Binary Model Evaluation

Validation Accuracy - **0.7181** Test Accuracy - **0.7289**



Multiple Binary-Classifiers Method

Pituitary Binary Classifier



Model Architecture

Same structure with Glioma Binary Model

3 – Convolutional Layer

3 – MaxPooling Layer

1 - Fully Connected Layer

Output Layer with 2 classes

Model Compiling

Model Training



Same paramaters with Multiclass CNN Model

.....


Binary-Classifiers Method

Pituitary Binary Model Evaluation

Validation Accuracy - **0.8275** Test Accuracy - **0.8675**



Ensemble Method



```

1 path='C:/Users/laptop/Desktop/BRAIN_TUMOR_STUDY/brain-tumor-classification-mri/Training/MultiClass/meningioma_tumor/m2 (15
2 img_process(path)
  
```

	multi_glioma	multi_meningioma	multi_no	multi_pituitary	model_meningioma	model_pituitary	model_glioma	Highest predicted class	True Class
0	0.043288	0.956708	4.925593e-08	0.000008	0.890799	0.993142	0.848488	pituitary	meningioma



05

MODEL EVALUATION

Model Evaluation

Evaluation Metrics



RECALL

Capturing Tumor Highly Important For Treatment

because...

PRECISION


Avoid improper treatment to a patient with no disease or different type of tumor

$$\text{Recall} = \frac{\text{correctly detected tumor (TP)}}{\text{all actual tumor (TP + FN)}}$$


$$\text{Precision} = \frac{\text{correctly detected tumor (TP)}}{\text{all tumor predicted (TP + FP)}}$$

Confusion Matrices

MULTICLASS MODEL




MENINGIOMA BINARY MODEL




Confusion Matrices

PITUITARY BINARY MODEL




GLIOMA BINARY MODEL



Recall and Precision Comparisons


MODEL EVALUATION

Meningioma Tumor Precision and Recall Metrics




Model_Multiclass Model_Meningioma Ensemble Models

Pituitary Tumor Precision and Recall Metrics




Model_Multiclass Model_Pituitary Ensemble Model

Glioma Tumor Precision and Recall Metrics



Model_Multiclass Model_Glioma Ensemble Models

Validation and Test Accuracy of Models





06

CONCLUSION

CONCLUSION

The highest Prediction and Recall of each tumor:

Meningioma Tumor Classification: **Recall:** Model_Meningioma - 0.91 **Precision:** Multiclass Model - 0.84

Pituitary Tumor Classification: **Recall:** Multiclass Model – 0.97 **Precision:** Model_Pituitary- 0.97

Glioma Tumor Classification: **Recall:** Multiclass Model – 0.82 **Precision:** Multiclass Model - 0.92

Almost all models with 90% precision and recall scores.

Multiclass CNN model with the highest either recall or precision score.

The highest Test Accuracy:

Multiclass CNN model with Test Accuracy 0.8815.

The Multiclass Model has the highest precision, recall and test accuracy.



07

RECOMMENDATION

RECOMMENDATION

Since Multiclass Model results almost 90% accuracy in detecting and classifying brain tumor we would recommend to use Multiclass Model.



08

FUTURE WORK

FUTURE WORK

To build more robust deep learning model by

- using more images
- applying GridSearch for hyperparameter optimization
- using GPUs to reduce training time